




Nitric 20%

Instructions



Initial Acid Passivation

1. After cleaning vessel, make a solution of 5.0 oz. of Nitric 20% to 1 gallon of water.
2. Circulate for 20-25 minutes at a temperature of 120°F
3. Drain vessel and **ALLOW TO AIR DRY** for 24 hours
4. Vessel must air dry in order to achieve successful passivation
5. Sanitize before use



Re-passivation

1. After cleaning vessel, make a solution of 3oz. of (Nitric 20% ,CIP CID #1832, CIP CID #1834) to 1 gallon of water
2. Circulate for 30 minutes at a temperature of 120°F
3. Drain vessel and **ALLOW TO AIR DRY**
4. Vessel must air dry in order to achieve successful passivation
5. Sanitize before use

Sales@stone-leigh.com
Stone-Leighbrewcleaners.com
@stoneleighbrewcleaner
508-465-0230

